


Eldor Wound Catheter Kit

Fenestrated catheters for local anaesthetic infusion

- Compatible with elastomeric and electronic pumps
- 16 gauge Tuohy needle - 80mm or 200mm
- Catheter with 10 lateral fenestrations
- High quality design and manufacture
- Cost effective



MODEL	DESCRIPTION
CSEN194	10cm fenestrations, 8cm introducer and luer lock connector
CSEN194.2	10cm fenestrations, 20cm introducer and luer lock connector

Clinical Applications

- Surgical Site Infusions
- Nerve Block Infusions




Customer Service

Australia T 1300 550 310 E info.au@admedus.com
New Zealand T 0800 226 926 E info.nz@admedus.com

Eldor Ultrasound Kit

Specialised kit for ultrasound guided nerve blocks



- Huber needle with echogenic properties
- Curved tip to facilitate ease of penetration, reducing risk of damage
- Two introducer lengths - 60 and 90 mm

Clinical Applications

- Continuous Peripheral Nerve Blocks
- Single Shot Nerve Block

MODEL	DESCRIPTION
CSEN 242	Tuohy Needle 18G x 90mm, Nerve Block Catheter, 3 lateral holes, Adapter, Filter, Thread Guide
CSEN 244.1	Tuohy Needle 18G x 60mm, Nerve Block Catheter, 3 lateral holes, Adapter, Filter, Thread Guide